

Human Bone Alkaline Phosphatase ELISA Kit | BALP elisa kit

Human Bone Alkaline Phosphatase ELISA Kit

CATALOG #
AAA10011

REACTIVITY
Human


PRODUCT SYNONYMS

Bone Alkaline Phosphatase; Human Bone Alkaline Phosphatase ELISA Kit; BALP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Quantitative Sandwich
Detection Range	1.56ng/ml-50ng/ml
Sensitivity	0.1ng/ml
Intra-assay Precision	Intra-assay CV (%) is less than 15%.
Inter-assay Precision	Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean ×100]. All CV% should be compared by concentration, not compared by OD values.
Preparation and Storage	Store all reagents at 2-8 degree C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BALP elisa kit

Background/Introduction: This Quantitative Sandwich ELISA kit is for lab reagent/research use only, not for drug, household, therapeutic or test applications! This kit is intended to be used for determine the level of BALP (hereafter termed "analyte") in undiluted original Human serum, plasma or tissue homogenates samples. For other sample types please contact tech support to determine compatibility with this assay. This kit is not suitable for assaying non-biological sources of substances.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.