

Human TIE2 ELISA Kit | TEK elisa kit

Human TIE2 ELISA Kit

CATALOG # AAA105548	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

TIE2; Human TIE2 ELISA Kit; VMCM; VMCM1; Angiopoietin-1 receptor; Endothelial tyrosine kinase; Tunica interna endothelial cell kinase; Tyrosine kinase with Ig and EGF homology domains-2; Tyrosine-protein kinase receptor TEK; Tyrosine-protein kinase receptor TIE-2; hTIE2; p140 TEK; CD202b; TEK elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	1124
Samples	Cell Culture Supernatants, Serum, Plasma
Assay Type	Sandwich
Detection Range	156 pg/ml - 10000 pg/ml
Sensitivity	5 pg/ml
Applicable Applications for TEK elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TEK elisa kit	Sandwich ELISA for quantitative detection of Human TIE2 concentrations.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	88758596
NCBI GeneID	7010

NCBI Accession #	NP_000450.2
NCBI GenBank Nucleotide #	NM_000459.4
UniProt Accession #	Q02763 ; Q5TCU2 ; Q8IV34 ; A8K6W0 ; B4DH20 ; B4DHD3 ; D3DRK5 ; E7EWI2
Molecular Weight	109,141 Da
NCBI Official Full Name	angiopoietin-1 receptor isoform 1
NCBI Official Synonym Full Names	TEK tyrosine kinase, endothelial
NCBI Official Symbol	TEK
NCBI Official Synonym Symbols	TIE2; VMCM; TIE-2; VMCM1; CD202B
NCBI Protein Information	angiopoietin-1 receptor; endothelial tyrosine kinase; tyrosine-protein kinase receptor TEK; tunica interna endothelial cell kinase; tyrosine-protein kinase receptor TIE-2; tyrosine kinase with Ig and EGF homology domains-2
UniProt Protein Name	Angiopoietin-1 receptor
UniProt Gene Name	TEK
UniProt Synonym Gene Names	TIE2; VMCM; VMCM1; hTIE2
UniProt Entry Name	TIE2_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/tek-human-elisa-kits/tie2/105548>