

Human CD80 ELISA Kit | CD80 elisa kit

Human CD80 ELISA Kit

CATALOG # AAA105673	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

CD80; Human CD80 ELISA Kit; CD28LG; CD28LG1; LAB7; T-lymphocyte activation antigen CD80; Activation B7-1 antigen; BB1; CTLA-4 counter-receptor B7.1; B7; CD80 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	288
Samples	Cell Culture Supernatants, Serum
Assay Type	Sandwich
Detection Range	62.5 pg/ml - 4000 pg/ml
Sensitivity	10 pg/ml
Applicable Applications for CD80 elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CD80 elisa kit	Sandwich ELISA for quantitative detection of Human CD80 concentrations.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	4885123
NCBI GeneID	941

NCBI Accession #	NP_005182.1
NCBI GenBank Nucleotide #	NM_005191.3
UniProt Accession #	P33681; Q5DTA9; Q5DTB0
Molecular Weight	18,761 Da
NCBI Official Full Name	T-lymphocyte activation antigen CD80
NCBI Official Synonym Full Names	CD80 molecule
NCBI Official Symbol	CD80
NCBI Official Synonym Symbols	B7; BB1; B7-1; B7.1; LAB7; CD28LG; CD28LG1
NCBI Protein Information	T-lymphocyte activation antigen CD80; activation B7-1 antigen; costimulatory factor CD80; CTLA-4 counter-receptor B7.1; B-lymphocyte activation antigen B7; costimulatory molecule variant IgV-CD80; CD80 antigen (CD28 antigen ligand 1, B7-1 antigen)
UniProt Protein Name	T-lymphocyte activation antigen CD80
UniProt Gene Name	CD80
UniProt Synonym Gene Names	CD28LG; CD28LG1; LAB7; B7
UniProt Entry Name	CD80_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/cd80/105673>