

Human REG4 ELISA Kit | REG4 elisa kit

Human REG4 ELISA Kit

CATALOG # AAA105700	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

REG4; Human REG4 ELISA Kit; GISP; RELP; Regenerating islet-derived protein 4; REG-4; Gastrointestinal secretory protein; REG-like protein; Regenerating islet-derived protein IV; Reg IV; REG4 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	158
Samples	Cell Culture Supernatants, Serum, Plasma
Assay Type	Sandwich
Detection Range	156 pg/ml - 10000 pg/ml
Sensitivity	10 pg/ml
Applicable Applications for REG4 elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for REG4 elisa kit	Sandwich ELISA for quantitative detection of Human REG4 concentrations.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	226823233
NCBI GeneID	83998

NCBI Accession #	NP_001152824.1
NCBI GenBank Nucleotide #	NM_001159352.1
UniProt Accession #	Q9BYZ8; Q8NER6; Q8NER7
Molecular Weight	14,993 Da
NCBI Official Full Name	regenerating islet-derived protein 4 isoform 1
NCBI Official Synonym Full Names	regenerating islet-derived family, member 4
NCBI Official Symbol	REG4
NCBI Official Synonym Symbols	GISP; RELP; REG-IV
NCBI Protein Information	regenerating islet-derived protein 4; REG-4; reg IV; REG-like protein; regenerating gene type IV; gastrointestinal secretory protein; regenerating islet-derived protein IV
UniProt Protein Name	Regenerating islet-derived protein 4
UniProt Gene Name	REG4
UniProt Synonym Gene Names	GISP; RELP; REG-4; Reg IV
UniProt Entry Name	REG4_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/reg4/105700>