

Mouse CXCL10 ELISA Kit | CXCL10 elisa kit

Mouse CXCL10 ELISA Kit

CATALOG # AAA105834	REACTIVITY Mouse	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

CXCL10; Mouse CXCL10 ELISA Kit; INP10; SCYB10; C-X-C motif chemokine 10; 10 kDa interferon gamma-induced protein; Gamma-IP10; IP-10; Small-inducible cytokine B10; CXCL10 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	98
Samples	Cell Culture Supernatants, Serum, Plasma
Assay Type	Sandwich
Detection Range	31.2 pg/ml - 2000 pg/ml
Sensitivity	1 pg/ml
Applicable Applications for CXCL10 elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CXCL10 elisa kit	Sandwich ELISA for quantitative detection of Mouse CXCL10 concentrations.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	10946576
NCBI GeneID	15945

NCBI Accession #	NP_067249.1
NCBI GenBank Nucleotide #	NM_021274.2
UniProt Accession #	P17515
Molecular Weight	10,789 Da
NCBI Official Full Name	C-X-C motif chemokine 10
NCBI Official Synonym Full Names	chemokine (C-X-C motif) ligand 10
NCBI Official Symbol	Cxcl10
NCBI Official Synonym Symbols	C7; IP10; CRG-2; INP10; IP-10; Ifi10; mob-1; Scyb10; gIP-10
NCBI Protein Information	C-X-C motif chemokine 10; gamma-IP10; interferon activated gene 10; small-inducible cytokine B10; interferon-gamma induced protein CRG-2; interferon-gamma-induced protein CRG-2; 10 kDa interferon gamma-induced protein; small inducible cytokine B subfamily (Cys-X-Cys), member 10
UniProt Protein Name	C-X-C motif chemokine 10
UniProt Gene Name	Cxcl10
UniProt Synonym Gene Names	Crg2; Ifi10; Inp10; Scyb10; Gamma-IP10; IP-10
UniProt Entry Name	CXL10_MOUSE

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/mouse-elisa-kits/cxcl10/105834>