

Human VEGFR3 ELISA Kit | FLT4 elisa kit

Human VEGFR3 ELISA Kit

CATALOG # AAA106081	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

VEGFR3; Human VEGFR3 ELISA Kit; Vascular endothelial growth factor receptor 3; VEGFR-3; Fms-like tyrosine kinase 4; FLT-4; Tyrosine-protein kinase receptor FLT4; FLT4 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	1298
Samples	Cell Culture Supernatants, Serum, Plasma
Assay Type	Sandwich
Detection Range	156 pg/ml - 10000 pg/ml
Sensitivity	10 pg/ml
Applicable Applications for FLT4 elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FLT4 elisa kit	Sandwich ELISA for quantitative detection of Human VEGFR3 concentrations.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	103472027
NCBI GeneID	2324

NCBI Accession #	NP_002011.2
NCBI GenBank Nucleotide #	NM_002020.4
UniProt Accession #	P35916 ; Q16067 ; Q86W07 ; Q86W08 ; A8K6L4 ; B5A926 ; AAB23636
Molecular Weight	93,174 Da
NCBI Official Full Name	vascular endothelial growth factor receptor 3 isoform 2
NCBI Official Synonym Full Names	fms-related tyrosine kinase 4
NCBI Official Symbol	FLT4
NCBI Official Synonym Symbols	PCL; FLT41; LMPH1A; VEGFR3
NCBI Protein Information	vascular endothelial growth factor receptor 3; FLT-4; VEGFR-3; soluble VEGFR3 variant 1; soluble VEGFR3 variant 2; soluble VEGFR3 variant 3; fms-like tyrosine kinase 4; tyrosine-protein kinase receptor FLT4
UniProt Protein Name	Vascular endothelial growth factor receptor 3
UniProt Gene Name	FLT4
UniProt Synonym Gene Names	VEGFR3; VEGFR-3; FLT-4
UniProt Entry Name	VGFR3_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/flt4-human-elisa-kits/vegfr3/106081>