

Solute carrier family 2, facilitated glucose transporter member 1 (SLC2A1) Recombinant Protein | SLC2A1 recombinant protein

Recombinant Human Solute carrier family 2, facilitated glucose transporter member 1 (SLC2A1), partial

CATALOG #
AAA113364

HOST
E Coli or Yeast or Baculovirus or Mammalian Cell

PRODUCT SYNONYMS

Solute carrier family 2; facilitated glucose transporter member 1 (SLC2A1); Recombinant Human Solute carrier family 2; partial; SLC2A1 recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	207-271aa; Partial
Sequence	CPESPRFLLI NRNEENRAKS VLKKLRGTAD VTHDLQEMKE ESRQMMREKK VTILELFRSP AYRQP
Species	Homo sapiens (Human)
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

NCBI AND UNIPROT INFORMATION

NCBI GI #	166795299
NCBI GeneID	6513
NCBI Accession #	NP_006507.2
NCBI GenBank Nucleotide #	NM_006516.2
UniProt Accession #	P11166; O75535; Q147X2; A8K9S6; B2R620; D3DPX0
Molecular Weight	15.2 kDa
NCBI Official Full Name	solute carrier family 2, facilitated glucose transporter member 1
NCBI Official Synonym Full Names	solute carrier family 2 member 1
NCBI Official Symbol	SLC2A1
NCBI Official Synonym Symbols	CSE; PED; DYT9; GLUT; DYT17; DYT18; EIG12; GLUT1; HTLVR; GLUT-1; SDCHCN; GLUT1DS
NCBI Protein Information	solute carrier family 2, facilitated glucose transporter member 1
UniProt Protein Name	Solute carrier family 2, facilitated glucose transporter member 1
UniProt Gene Name	SLC2A1
UniProt Synonym Gene Names	GLUT1; GLUT-1

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

