

Outer capsid protein sigma-1 (S1) Recombinant Protein | S1 recombinant protein

Recombinant Reovirus type 3 Outer capsid protein sigma-1 (S1)

CATALOG # AAA114323	HOST E Coli or Yeast or Baculovirus or Mammalian Cell
-------------------------------	---


PRODUCT SYNONYMS

Outer capsid protein sigma-1 (S1); Recombinant Reovirus type 3 Outer capsid protein sigma-1 (S1); Outer capsid protein sigma-1; Sigma1; Cell attachment protein; Hemagglutinin; S1 recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	1-455. Full Length
Sequence	MDPRLREEVV RLIALTSDN GASLSKGLES RVSALEKTSQ IHSDTILRIT QGLDDANKRI IALEQSRDDL VASVSDAQLA ISRLESSIGA LQTVVNGLDS SVTQLGARVG QLETGLAELR VDHDNLVARV DTAERNIGSL TTELSTLTLR VTSIQADFES RISTLERTAV TSAGAPLSIR NNRMTMGLND GLTSGNNLA IRLPGNTGLN IQNGGLQFRF NTDQFQIVNN NLTKTTVFD SINSRIGATE QSYVASAVTP LRLNSSTKVL DMLIDSSTLE INSSGQLTVR STSPNLRYPI ADVSGGIGMS PNYRFRQSMW IGIVSYSGSG LNWRVQVNSD IFIVDDYIHI CLPAFDGFSI ADGGDLSLNF VTGLLPPLLT GDTEPAFHND VVTYGAQTVA IGLSSGGAPQ YMSKNLWVEQ WQDGVRLRV EGGGSITHSN SKWPAMTVSY PRSFT
Species	Reovirus type 3 (strain Dearing) (T3D) (Mammalian orthoreovirus 3)
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

IMAGES



Application Data

NCBI AND UNIPROT INFORMATION

NCBI GI #	205830900
UniProt Accession #	P03528; Q85668; A4ZY26
Molecular Weight	49,095 Da
NCBI Official Full Name	Outer capsid protein sigma-1
UniProt Protein Name	Outer capsid protein sigma-1
UniProt Gene Name	S1

UniProt Synonym Gene Names	Sigma1
UniProt Entry Name	SIGM1_REOVD

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/s1-recombinant-protein-e-coli/114323>