

Probable periplasmic serine endoprotease DegP-like (htrA) Recombinant Protein | htrA recombinant protein

Recombinant Bartonella henselae Probable periplasmic serine endoprotease DegP-like (htrA)

CATALOG #
AAA114425

HOST
E Coli or Yeast or Baculovirus or Mammalian Cell

PRODUCT SYNONYMS

Probable periplasmic serine endoprotease DegP-like (htrA); Recombinant Bartonella henselae Probable periplasmic serine endoprotease DegP-like (htrA); htrA recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	19-503aa; Full Length of Mature Protein
Sequence	LETALFFSGC GSSLWTTKAH ANSVFSSLMQ QQGFADIVSQ VKPAVVSQV KSNKKKKEWF FSDFFSTPGF DQLPDQHPLK KFFQDFYNRD KPSNKSQRS HRLRPIAFGS GFFISSDGYI VTNNHVISDG TSYAVVLDDG TELNAKLIGT DPRTDLAVLK VNEKRKFSYV DFGDDSKLRV GDWVVAIGNP FGLGGTVTAG IVSARGRDIG TGVYDDFIQI DAAVNRGNSG GPTFDLNGKV VGVNTAIFSP SGGNVGIAFA IPAATAKQVV QQLIEKGLVQ RGWLGVQIQP VTKEISDSIG LKEAKGALIT DPLKGPAKA GIKAGDVIIS VNGEKINDVR DLAKRIANMS PGETVTLGWV KSGKEENIKV KLSMPEDEN MKDGSKYSNE HGNSDETLED YGLIVAPSDD GVGLVVTDVD PDSDAADKGI RPDGVIVTVN NKSVKKVS DI TDTIKNAQKL GRKAILLQVR TNDQNRFVAL PIFKK
Species	Bartonella henselae (strain 49882 / Houston 1) (Rochalimaea henselae)
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

NCBI AND UNIPROT INFORMATION

NCBI GI #	499493768
NCBI Accession #	WP_011180408.1
NCBI GenBank Nucleotide #	WP_011180408.1
UniProt Accession #	P54925; Q6G491
Molecular Weight	54,113 Da
NCBI Official Full Name	DegP-like serine endoprotease
UniProt Protein Name	Probable periplasmic serine endoprotease DegP-like
UniProt Gene Name	htrA

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.