

Outer capsid protein lambda-2 (L2) Recombinant Protein | L2 recombinant protein

Recombinant Reovirus type 3 Outer capsid protein lambda-2 (L2), partial

CATALOG #
AAA115084

HOST
E Coli or Yeast or Baculovirus or Mammalian Cell

PRODUCT SYNONYMS

Outer capsid protein lambda-2 (L2); Recombinant Reovirus type 3 Outer capsid protein lambda-2 (L2); partial; L2 recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	1-200. partial
Sequence	MANVWGVRLA DSLSSPTIET RTRQYTLHDL CSDLDPANPGR EPWKPLRNQR TNNIVAVQLF RPLQGLVLDL QLYGFPGAFD DWERFMREKL RVLKYEVLRI YPISNYSNEH VNVFVANALV GAFLSNQAFY DLLPLLIIND TMIGDLLGTG ASLSQFFQSH GDVLEVAAGR KYLQMENYSN DDDDPPLFAK DLSDYAKAFY
Species	Reovirus type 3 (strain Dearing) (T3D) (Mammalian orthoreovirus 3)
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

NCBI AND UNIPROT INFORMATION

NCBI GI #	205737936
UniProt Accession #	P11079; A4ZY21
Molecular Weight	143,976 Da
NCBI Official Full Name	Outer capsid protein lambda-2
UniProt Protein Name	Outer capsid protein lambda-2
UniProt Gene Name	L2
UniProt Synonym Gene Names	Lambda2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/l2-recombinant-protein-e-coli/115084>