

Fusion glycoprotein F0 (F) Recombinant Protein | F recombinant protein

Recombinant Hendra virus Fusion glycoprotein F0 (F), partial

CATALOG # AAA115503	HOST E Coli or Yeast or Baculovirus or Mammalian Cell
-------------------------------	---

PRODUCT SYNONYMS

Fusion glycoprotein F0 (F); Recombinant Hendra virus Fusion glycoprotein F0 (F); partial; F recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	27-494. Partial, provide the complete extracellular domain.
Sequence	ILHYEKLSKI GLVKGITRKY KIKSNPLTKD IVIKMIPNVS NVSKCTGTVM ENYKSRLTGI LSPIKGAIEL YNNNTHDLVG DVKLAGVVMA GIAIGIATAA QITAGVALYE AMKNADNINK LKSSIESTNE AVVKLQETAE KTVYVLTALQ DYINTNLVPT IDQISCKQTE LALDLALSKY LSDLLFVFGP NLQDPVSNSM TIQAISQAFG GNYETLLRTL GYATEDFDDL LESDSIAGQI VYVDLSSYYI IVRVYFPILT EIQQAYVQEL LPVSFNNDNS EWISIVPNFV LIRNTLISNI EVKYCLITKK SVICNQDYAT PMTASVRECL TGSTDKCPRE LVVSSHVPRF ALSGGVLFAN CISVTCQCQT TGRAISQSGE QTLLMIDNTT CTTVVLGNII ISLGKYLGSI NYNSESIAVG PPVYTDKVDI SSQISSMNQS LQQSKDYIKE AQKILDTVNP SLISMLSM
Sequence Length	494
Species	Hendra virus (isolate Horse/Australia/Hendra/1994)
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

NCBI AND UNIPROT INFORMATION

NCBI GI #	29468606
NCBI GeneID	1446467
NCBI Accession #	NP_047111.2
NCBI GenBank Nucleotide #	NC_001906.3
UniProt Accession #	O89342 ; Q66761
Molecular Weight	59,795 Da
NCBI Official Full Name	fusion
NCBI Official Symbol	F
NCBI Protein Information	fusion
UniProt Protein Name	Fusion glycoprotein F0
UniProt Gene Name	F
UniProt Synonym Gene Names	Protein F

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/f-recombinant-protein-e-coli/115503>