

Gamma-aminobutyric acid receptor subunit pi (GABRP) Recombinant Protein | GABRP recombinant protein

Recombinant Human Gamma-aminobutyric acid receptor subunit pi (GABRP), partial

CATALOG #
AAA115845

HOST
E Coli or Yeast or Baculovirus or Mammalian Cell

PRODUCT SYNONYMS

Gamma-aminobutyric acid receptor subunit pi (GABRP); Recombinant Human Gamma-aminobutyric acid receptor subunit pi (GABRP); partial; GABRP recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	17-242aa; partial
Sequence	ERMCIQGSQF NVEVGRSDKL SLPGFENLTA GYNKFLRPNF GGEPVQIALT LDIASISSIS ESNMDYTATI YLRQRWMDQR LVFEGNKSFT LDARLVEFLW VPDYIYESK KSFLHEVTVG NRLIRLFSNG TVLYALRITT TVACNMDLSK YPMDTQTCKL QLESWGYDGN DVEFTWLRGN DSVRGLEHLR LAQYTIERYF TLVTRSQQET GNYTRLVLQF ELRRNV
Species	Homo sapiens (Human)
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

NCBI AND UNIPROT INFORMATION

NCBI GI #	635172858
NCBI GeneID	2568
NCBI Accession #	NP_001278914.1
NCBI GenBank Nucleotide #	NM_001291985.1
UniProt Accession #	O00591 ; Q32MJ1 ; A8KA36 ; D3DQL2
Molecular Weight	50,640 Da
NCBI Official Full Name	gamma-aminobutyric acid receptor subunit pi isoform 2
NCBI Official Synonym Full Names	gamma-aminobutyric acid type A receptor pi subunit
NCBI Official Symbol	GABRP
NCBI Protein Information	gamma-aminobutyric acid receptor subunit pi
UniProt Protein Name	Gamma-aminobutyric acid receptor subunit pi
UniProt Gene Name	GABRP

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

