

Methyl-accepting chemotaxis protein I (tsr) Recombinant Protein | tsr recombinant protein

Recombinant Escherichia coli Methyl-accepting chemotaxis protein I (tsr), partial

CATALOG # AAA115931	HOST E Coli or Yeast or Baculovirus or Mammalian Cell
-------------------------------	---


PRODUCT SYNONYMS

Methyl-accepting chemotaxis protein I (tsr); Recombinant Escherichia coli Methyl-accepting chemotaxis protein I (tsr); partial; tsr recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	211-551aa; Partial
Sequence	WFGIKASLVA PMNRLIDSIR HIAGGDLVKP IEVDGSNEMG QLAESLRHMQ GELMRTVGDV RNGANAIYSG ASEIATGNND LSSRTEQQAA SLEETAASME QLTATVKQNA ENARQASHLA LSASETAQRG GKVDNVDVQT MRDISTSSQK IADIISVIDG IAFQTNILAL NAAVEAARAG EQGRGFAVVA GEVRNLAQRS AQAAREIKSL IEDSVGKVDV GSTLVESAGE TMAEIVSAVT RVTDIMGEIA SASDEQSRGI DQVGLAVAEM DRVTQQNAAL VEESAAAAAA LEEQASRLTE AVAVFRIQQQ QRETSAVVKT VTPAAPRKMA VADSEENWET F
Species	Escherichia coli (strain K12)
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

IMAGES



SDS-PAGE

NCBI AND UNIPROT INFORMATION

NCBI GI #	16132176
NCBI GeneID	948884
NCBI Accession #	NP_418775.1
NCBI GenBank Nucleotide #	NC_000913.3

UniProt Accession #	P02942; P76817; Q2M5W1
Molecular Weight	38 kDa
NCBI Official Full Name	methyl-accepting chemotaxis protein I, serine sensor receptor
NCBI Official Symbol	tsr
NCBI Official Synonym Symbols	cheD; ECK4345; JW4318
NCBI Protein Information	methyl-accepting chemotaxis protein I, serine sensor receptor
UniProt Protein Name	Methyl-accepting chemotaxis protein I
UniProt Gene Name	tsr
UniProt Synonym Gene Names	cheD; MCP-I

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/tsr-recombinant-protein-e-coli/115931>