

P2RY11 cDNA clone

P2RY11 cDNA Clone

CATALOG #
AAA116328

PRODUCT SYNONYMS

P2RY11; P2RY11 cDNA Clone; P2RY11 cDNA clone

PRODUCT SPECIFICATIONS

Sequence	atggcagcca acgtctcggg tgccaagtcc tgccctgcc acttcttggc agctgccgac gacaaactca gtgggttcca gggggacttc ctgtggcca tactggtggt tgagttcctg gtggccgtgg ccagcaatgg cctggccctg taccgttca gcatccgaa gcagcgcca tggcaccg ccgtggtt ctctgtccag ctggcagtca gcgacctgt ctgcgtctg acgtgcccc cgctggcgc ctacctat cccccaagc actggcgta tggggaggcc gcgtgccgct tggagcgct cctttcacc tgcaacctgc tgggcagcgt catctcatc acctgcatc gcctcaaccg ctacctggc atcgtgcacc cttcttcgc ccgaagccac ctgcgaccca agcacgcctg ggccgtgagc gctgccggt gggctctggc cgccctgctg gccatgcca cactcagct ctccacctg aagaggcgc agcagggggc gggcaactgc agcgtggcca ggcccaggc ctgcatcaag tgtctggga cagcagacca cgggctggcg gcctacagag cgtatagcct ggtgtggcg ggggtggct gcggcctgcc gctgctgct acgtggcag cctacggcg cctcggggcg gccgtgtac gcagcccagg catgactgtg gccgagaagc tgcgtgtggc agcgttggg gccagtggg tggccctcta cgccagctcc tatgtgcct accacatcat cggggtgtc aacgtggatg ctggcgcg cgtagcacc cgctgccga gcttgcaga catagcccag gccacagcag cctggagct ggggcccac gtgggtacc aggtgatcg gggcctatg ccctggcct tctgttcca cctctact tacatggcg cagtcccag cctgggtgc tgcgtccgac actgcccgg ctacaggac agctggaacc cagaggacgc caagagcact ggccaagccc tgccctcaa tgccacagc gccctaaac cgtcagagcc ccagtccgt gagctgagcc aatga
Sequence Length	1125
Vector	Please Inquire

NCBI AND UNIPROT INFORMATION

NCBI GI #	49258118
NCBI GeneID	5032
Molecular Weight	40,345 Da
NCBI Official Full Name	Homo sapiens purinergic receptor P2Y, G-protein coupled, 11, mRNA
NCBI Official Synonym Full Names	purinergic receptor P2Y11
NCBI Official Symbol	P2RY11
NCBI Official Synonym Symbols	P2Y11
NCBI Protein Information	P2Y purinoceptor 11
UniProt Protein Name	P2Y purinoceptor 11
UniProt Gene Name	P2RY11
UniProt Synonym Gene Names	P2Y11
UniProt Entry Name	P2Y11_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.