

LPS-assembly protein lptD (lptD) Recombinant Protein | lptD recombinant protein

Recombinant Francisella tularensis subsp. tularensis LPS-assembly protein lptD (lptD), partial

CATALOG #
AAA117620

HOST
E Coli or Yeast or Baculovirus or Mammalian Cell

PRODUCT SYNONYMS

LPS-assembly protein lptD (lptD); Recombinant Francisella tularensis subsp. tularensis LPS-assembly protein lptD (lptD); partial; lptD recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	25-484aa; partial
Sequence	ARIMSNPIK EDWQCKVVDG EWSCKRAKKP KSVFDKCLK TEKEKALADD LAWVKKPSYF VGGYYSNDNQ FTKALCESKK TDLSYEKSEF DNYGTLIASG NVQVLQCDQE LYGNNAIINL NSNNSAIRSL VMAGDVIVKQ PSTGIVIRTT ELDADMNNGT YSTGEAYFRL AREMPKTRIY DKEHFSGYLR GYAKTFKKES SGDVLSDGY ITSGDPYDNA WKITGNINIDI DTNTHMAYVK NGYFEIQDIP VMYIPYFSHP IDRRRRSGFL YPGFVQNANS GIGISVPYYF NLAPNYDLML QSVIWSQRGI IENGTFRYMT KYFQQQFEFS LVPYDFKEGK MRGSFTLSTT GQYENINTNF KYEYVSDQNY YNDFSAGNVN LVTKLLDRE FDLTYTNDYV DSGLTVLDYG VVNPLLVDN TPYAKLPEVK LNLTS DGYP DYLTLSAQTL NTFFYKTAGP
Species	Francisella tularensis subsp. tularensis (strain SCHU S4 / Schu 4)
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

NCBI AND UNIPROT INFORMATION

NCBI GI #	489113349
NCBI Accession #	WP_003023202.1
NCBI GenBank Nucleotide #	WP_003023202.1
UniProt Accession #	Q5NHI7
Molecular Weight	98,401 Da
NCBI Official Full Name	LPS-assembly protein LptD
UniProt Protein Name	LPS-assembly protein LptD
UniProt Gene Name	lptD

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.