

Rabbit Peptide deformylase Polyclonal Antibody | anti-def antibody

Rabbit anti-Escherichia coli Peptide deformylase polyclonal Antibody


CATALOG # AAA118989	HOST Rabbit	REACTIVITY Escherichia coli	APPLICATIONS Western Blot, ELISA
-------------------------------	-----------------------	---------------------------------------	--

PRODUCT SYNONYMS

Peptide deformylase; Polyclonal Antibody; Rabbit anti-Escherichia coli Peptide deformylase polyclonal Antibody; Peptide deformylase (PDF) (EC 3.5.1.88) (Polypeptide deformylase); def; fms; anti-def antibody

PRODUCT SPECIFICATIONS

Host	Rabbit
Reactivity	Escherichia coli
Clonality	Polyclonal
Isotype	IgG
Purity/Purification	>95%, Protein G purified
Form/Format	Liquid Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Sequence Length	169
Applicable Applications for anti-def antibody	WB (Western Blot), ELISA
Conjugate	Non-conjugated
Immunogen	Recombinant Escherichia coli Peptide deformylase protein (2-169aa)
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	def
Preparation and Storage	Upon receipt, store at -20 degrees C or -80 degrees C. Avoid repeated freeze.



WB (Western Blot)
Western blotAll lanes: def antibody at 2ug/ml + DH5 alpha whole cell lysateSecondaryGoat polyclonal to rabbit at 1/10000 dilutionPredicted band size: 20kDaObserved band size: 20kDa

NCBI AND UNIPROT INFORMATION

NCBI GI #	16131166
NCBI GeneID	947780
NCBI Accession #	NP_417745.1
NCBI GenBank Nucleotide #	NC_000913.3
UniProt Accession #	P0A6K3; P27251; Q2M6V1; CAA54826
NCBI Official Full Name	peptide deformylase
NCBI Official Symbol	def
NCBI Official Synonym Symbols	ECK3273; fms; JW3248
NCBI Protein Information	peptide deformylase
UniProt Protein Name	Peptide deformylase
UniProt Gene Name	def
UniProt Synonym Gene Names	fms; PDF
UniProt Entry Name	DEF_ECOLI

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.