

Teprotumumab ELISA Kit

Teprotumumab ELISA Kit

CATALOG #
AAA120185

PRODUCT SYNONYMS

Teprotumumab; Teprotumumab ELISA Kit; RO4858696-000; Teprotumumab elisa kit

PRODUCT SPECIFICATIONS

Assay Type	Quantitative Competitive
Samples	Serum and plasma
Detection Range	62.5 ng/mL - 4,000 ng/mL
Intra-assay Precision	Intra-Assay Precision (Precision within an assay): <20%. Three samples of known concentration were tested sixteen times on one plate to assess intra-assay precision.
Inter-assay Precision	Inter-Assay Precision (Precision between assays): <20%. Three samples of known concentration were tested in twenty four separate assays to assess inter-assay precision.
Preparation and Storage	Store at 2-8 degree C The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.

ADDITIONAL INFORMATION

Related Product Information for Teprotumumab elisa kit	Principle of the Assay: This assay employs the quantitative competitive enzyme immunoassay technique. Recombinant human IGF1R antigen has been pre-coated onto a microplate. Standards or samples are premixed with biotin-labeled antibody and then pipetted into the wells. Teprotumumab in the sample competitively binds to the pre-coated protein with biotin-labeled Teprotumumab. After washing away any unbound substances, Streptavidin-HRP is added to the wells. Following a wash to remove any unbound enzyme reagent, a substrate solution is added to the wells and color develops in inversely proportion to the amount of Teprotumumab bound in the initial step. The color development is stopped and the intensity of the color is measured. Intended Uses: Used for the quantitative determination of Teprotumumab concentration in serum and plasma.
Product Categories/Family for Teprotumumab elisa kit	Assay Kits

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.comOnline datasheet: <https://www.aaabiotech.com/elisa-kits/teprotumumab/120185>