

Porcine tyrosine hydroxylase (TH) ELISA Kit | TH elisa kit

Porcine tyrosine hydroxylase (TH) ELISA Kit

CATALOG #
AAA121354

REACTIVITY
Porcine


PRODUCT SYNONYMS

tyrosine hydroxylase (TH); Porcine tyrosine hydroxylase (TH) ELISA Kit; TH elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Sequence Length	106
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	1ng/ml - 300ng/ml
Sensitivity	0.53ng/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TH elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Porcine Tyrosine Hydroxylase (also known as TH) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Porcine TH antibody. TH present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Porcine TH Antibody is added and binds to TH in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated TH antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Porcine TH. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

NCBI AND UNIPROT INFORMATION

NCBI GI #	986909
NCBI GeneID	7054
Molecular Weight	58,160 Da
NCBI Official Full Name	tyrosine hydroxylase, partial
NCBI Official Synonym Full Names	tyrosine hydroxylase
NCBI Official Symbol	TH
NCBI Official Synonym Symbols	TYH; DYT14; DYT5b
NCBI Protein Information	tyrosine 3-monooxygenase
UniProt Protein Name	Tyrosine 3-monooxygenase
UniProt Gene Name	TH
UniProt Synonym Gene Names	TYH; TH

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/porcine-elisa-kits/tyrosine-hydroxylase-th/121354>