

## Fish Cyclic adenosine monophosphate (cAMP) ELISA Kit | cAMP elisa kit

Fish cyclic adenosine monophosphate (cAMP) ELISA Kit

CATALOG #  
**AAA121846**

REACTIVITY  
**Fish**


### PRODUCT SYNONYMS

Cyclic adenosine monophosphate (cAMP); Fish cyclic adenosine monophosphate (cAMP) ELISA Kit; cAMP elisa kit

### PRODUCT SPECIFICATIONS

|                         |                                                                          |
|-------------------------|--------------------------------------------------------------------------|
| Reactivity              | Fish                                                                     |
| Samples                 | Serum, plasma, cell culture supernates, cell lysates, tissue homogenates |
| Assay Type              | Quantitative Sandwich                                                    |
| Detection Range         | 1nmol/L - 300nmol/L                                                      |
| Sensitivity             | 0.54nmol/L                                                               |
| Intra-Assay             | CV<8%                                                                    |
| Inter-Assay             | CV<10%                                                                   |
| Preparation and Storage | Store at 2 to 8 degree C for 6 months.                                   |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for cAMP elisa kit

**Intended Uses:** This sandwich kit is for the accurate quantitative detection of Fish Cyclic Adenosine Monophosphate (also known as cAMP) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

**Principles of the Assay:** This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Fish cAMP antibody. cAMP present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Fish cAMP Antibody is added and binds to cAMP in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated cAMP antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Fish cAMP. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)

Online datasheet: <https://www.aaabiotech.com/fish-elisa-kits/cyclic-adenosine-monophosphate-camp/121846>