

Canine Ferritin, FE ELISA Kit | FE elisa kit

Canine Ferritin, FE ELISA Kit

CATALOG #
AAA123883

REACTIVITY
Canine


PRODUCT SYNONYMS

Ferritin; FE; Canine Ferritin; FE ELISA Kit; FE elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Assay Type	Quantitative Sandwich
Samples	Serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids
Detection Range	1-400ng/ml
Sensitivity	0.59ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FE elisa kit

Principle of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Canine FE antibody. FE present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Canine FE Antibody is added and binds to FE in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated FE antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Canine FE. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

Intended Uses: This sandwich kit is for the accurate quantitative detection of Canine Ferritin (also known as FE) in serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.