

Mouse anti-Human Ig Lambda Light Chain Monoclonal Antibody | anti-IgI antibody

Anti-Hu Ig Lambda Light Chain PE, Immunoglobulin lambda, IgI

CATALOG # AAA128258	HOST Mouse	REACTIVITY Human	APPLICATIONS Flow Cytometry, Functional Assay
-------------------------------	----------------------	----------------------------	---

PRODUCT SYNONYMS

Ig Lambda Light Chain; Monoclonal Antibody; Anti-Hu Ig Lambda Light Chain PE; Immunoglobulin lambda; IgI; anti-IgI antibody

PRODUCT SPECIFICATIONS

Host	Mouse
Reactivity	Human
Clonality	Monoclonal
Isotype	Ig Lambda Light Chain, IgG1 kappa
Clone Number	1-155-2
Specificity	The mouse monoclonal antibody 1-155-2 recognizes lambda light chains (22.5 kDa) of human immunoglobulin.
Purity/Purification	The purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. The conjugate is purified by size-exclusion chromatography.
Form/Format	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Sequence Length	125
Applicable Applications for anti-IgI antibody	FCM/FACS (Flow Cytometry)
Immunogen	Human IgA1 lambda myeloma protein
Excitation laser	Blue (488nm)
Conjugation	PE
Preparation	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography
Preparation and Storage	Store in the dark at 2-8 degree. Do not freeze. Avoid prolonged exposure to light.

NCBI AND UNIPROT INFORMATION

NCBI GI #	1165318
NCBI Official Full Name	Ig lambda light chain, partial

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/monoclonal/ig-lambda-light-chain-antibody-clone-1-155-2-flow-cytometry-elisa/128258>