

Human Cleavage And Polyadenylation Specific Factor 1 (CPSF1) ELISA Kit | CPSF1 elisa kit

Human Cleavage And Polyadenylation Specific Factor 1 (CPSF1) ELISA Kit

CATALOG # AAA136349	REACTIVITY Human
-------------------------------	----------------------------


PRODUCT SYNONYMS

Cleavage And Polyadenylation Specific Factor 1 (CPSF1); Human Cleavage And Polyadenylation Specific Factor 1 (CPSF1) ELISA Kit; CPSF160; Cleavage and polyadenylation specificity factor 160 kDa subunit; CPSF1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Cleavage And Polyadenylation Specific Factor 1 (CPSF1). No significant cross-reactivity or interference between Cleavage And Polyadenylation Specific Factor 1 (CPSF1) and analogues was observed.
Sequence Length	377
Assay Type	Double-antibody Sandwich
Samples	Tissue homogenates and other biological fluids
Detection Range	0.312-20ng/mL
Sensitivity	0.129ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Cleavage And Polyadenylation Specific Factor 1 (CPSF1) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Cleavage And Polyadenylation Specific Factor 1 (CPSF1) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CPSF1 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Cleavage And Polyadenylation Specific Factor 1 (CPSF1). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Cleavage And Polyadenylation Specific Factor 1 (CPSF1). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Cleavage And Polyadenylation Specific Factor 1 (CPSF1), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Cleavage And Polyadenylation Specific Factor 1 (CPSF1) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for CPSF1 elisa kit	Metabolic pathway

NCBI AND UNIPROT INFORMATION

NCBI GI #	119602515
NCBI GeneID	29894
UniProt Accession #	Q10570; Q96AF0
Molecular Weight	160,884 Da
NCBI Official Full Name	cleavage and polyadenylation specific factor 1, 160kDa, isoform CRA_b
NCBI Official Synonym Full Names	cleavage and polyadenylation specific factor 1, 160kDa
NCBI Official Symbol	CPSF1
NCBI Official Synonym Symbols	CPSF160; P/cl.18; HSU37012
NCBI Protein Information	cleavage and polyadenylation specificity factor subunit 1
UniProt Protein Name	Cleavage and polyadenylation specificity factor subunit 1
UniProt Gene Name	CPSF1
UniProt Synonym Gene Names	CPSF160; CPSF 160 kDa subunit
UniProt Entry Name	CPSF1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.