

Bovine Fetuin B (FETUB) ELISA Kit | FETUB elisa kit

Bovine Fetuin B (FETUB) ELISA Kit

CATALOG #
AAA137609

REACTIVITY
Bovine


PRODUCT SYNONYMS

Fetuin B (FETUB); Bovine Fetuin B (FETUB) ELISA Kit; FETU-B; 16G2; Gugu; IRL685; Fetuin-like protein IRL685; FETUB elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Bovine
Specificity	This assay has high sensitivity and excellent specificity for detection of Fetuin B (FETUB). No significant cross-reactivity or interference between Fetuin B (FETUB) and analogues was observed.
Sequence Length	382
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma and other biological fluids
Detection Range	31.25-2000ng/mL
Sensitivity	11.7ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Fetuin B (FETUB) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Fetuin B (FETUB) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FETUB elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Fetuin B (FETUB). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Fetuin B (FETUB). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Fetuin B (FETUB), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Fetuin B (FETUB) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for FETUB elisa kit	Metabolic pathway; Infection immunity; Endocrinology; Bone metabolism

NCBI AND UNIPROT INFORMATION

NCBI GI #	58331240
NCBI GeneID	26998
NCBI Accession #	NP_055190.2
NCBI GenBank Nucleotide #	NM_014375.2
UniProt Accession #	Q9UGM5; Q1RMZ0; Q5J876; Q6DK58; Q6GRB6; Q9Y6Z0; B2RCW6; E9PG06
Molecular Weight	37,823 Da
NCBI Official Full Name	fetuin-B isoform 1
NCBI Official Synonym Full Names	fetuin B
NCBI Official Symbol	FETUB
NCBI Official Synonym Symbols	16G2; Gugu; IRL685
NCBI Protein Information	fetuin-B
UniProt Protein Name	Fetuin-B
UniProt Gene Name	FETUB
UniProt Entry Name	FETUB_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/bovine-ovine-elisa-kits/fetuin-b-fetub/137609>