

Mouse Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9) ELISA Kit | PCSK9 elisa kit

Mouse Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9) ELISA Kit

CATALOG #
AAA138407

REACTIVITY
Mouse


PRODUCT SYNONYMS

Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9); Mouse Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9) ELISA Kit; FH3; HCHOLA3; NARC1; Hypercholesterolemia; Autosomal Dominant 3; Neural apoptosis-regulated convertase 1; Proprotein convertase 9; Subtilisin/kexin-like protease PC9; PCSK9 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9). No significant cross-reactivity or interference between Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9) and analogues was observed.
Sequence Length	196
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma, Tissue homogenates and other biological fluids
Detection Range	15.6-1000pg/mL
Sensitivity	6.0pg/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PCSK9 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for PCSK9 elisa kit	Enzyme & Kinase

NCBI AND UNIPROT INFORMATION

NCBI GI #	31317307
NCBI GeneID	255738
NCBI Accession #	NP_777596.2
NCBI GenBank Nucleotide #	NM_174936.3
UniProt Accession #	Q8NBP7; Q5PSM5; Q5SZQ2; A8T640; C0JYY9
Molecular Weight	20,827 Da
NCBI Official Full Name	proprotein convertase subtilisin/kexin type 9 preproprotein
NCBI Official Synonym Full Names	proprotein convertase subtilisin/kexin type 9
NCBI Official Symbol	PCSK9
NCBI Official Synonym Symbols	FH3; PC9; NARC1; LDLCQ1; NARC-1; HCHOLA3
NCBI Protein Information	proprotein convertase subtilisin/kexin type 9
UniProt Protein Name	Proprotein convertase subtilisin/kexin type 9
UniProt Gene Name	PCSK9
UniProt Synonym Gene Names	NARC1; NARC-1; PC9
UniProt Entry Name	PCSK9_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.