

Human Immunoglobulin G3 (IgG3) DIY ELISA Kit | IgG3 diy elisa kit

Human Immunoglobulin G3 (IgG3) ELISA Kit DIY Materials

CATALOG #
AAA147650

REACTIVITY
Human


PRODUCT SYNONYMS

Immunoglobulin G3 (IgG3); Human Immunoglobulin G3 (IgG3) ELISA Kit DIY Materials; IGHG3; HDC; Heavy chain disease protein; Ig Gamma-3 Chain C Region; Immunoglobulin Heavy Constant Gamma 3; G3m Marker; Immunoglobulin Gm3; IgG3 diy elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	The Antibody Pairs in the kit have high sensitivity and excellent specificity for detection of IFN γ .
Sequence Length	377
Assay Type	Sandwich
Detection Range	3.12-200pg/mL.
Preparation and Storage	Antibodies, Standard and Streptavidin-HRP should be stored at -20 degree C. TMB should be stored at 4 degree C. 96-well Plate could be stored at room temperature. The contents are valid for twelve months. They are stable for one month after opening when stored at 4 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for IgG3 diy elisa kit

Intended Uses: This 'Do it Yourself (DIY)' assay kit contains materials for developing your own immunoassays. This may include the development of sandwich ELISA to measure IFN γ in gallus serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids in vitro. This kit contains sufficient materials for preparation of at least five 96-well plates. Sample Preparation, Assay Protocol and Detection Range for Reference are provided as suggestions only. Researchers should optimize the use of kit materials and protocols for their own model system and determine if the kit is suitable for their samples and immunoassay.

Product Categories/Family for IgG3 diy elisa kit

Infection immunity; Immune molecule

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.