

Human 5'-Nucleotidase, Cytosolic III (NT5C3) Antibody Pair Kit | NT5C3 antibody pair kit

Human 5'-Nucleotidase, Cytosolic III (NT5C3) Antibody Pair Kit (with Standard)

CATALOG # AAA148907	REACTIVITY Human	APPLICATIONS ELISA, ELISA
-------------------------------	----------------------------	-------------------------------------

PRODUCT SYNONYMS

5'-Nucleotidase; Cytosolic III (NT5C3); Human 5'-Nucleotidase; Cytosolic III (NT5C3) Antibody Pair Kit (with Standard); UMPH1; PSN1; PN-I; UMPH; P5'N-1; p36; cN-III; Pyrimidine 5'-nucleotidase 1; Uridine 5'-monophosphate hydrolase 1; NT5C3 antibody pair kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, Plasma, Tissue Homogenates, Cell Lysates, Cell Culture Supernates and Other Biological Fluids
Detection Range	31.2-2000ng/mL
Applicable Applications for NT5C3 antibody pair kit	Luminex, ELISPOT, CLIA (Chemiluminescence Immunoassay), ELISA
Standard	Recombinant NT5C3 (Lyophilized, 1 vial-10ug)
Capture Antibody	Rabbit polyclonal antibody (Liquid, 1 vial-400mug/800muL)
Detection Antibody	Biotinylated Rabbit polyclonal antibody (Liquid, 1 vial-100mug/500muL)
Producte Usage	Raw materials can be used for preparation of 5x96T or 10x96T ELISA Kits.
Antibody Support Pack	Antibody Support Pack (Buffer Pack) is sold separately as

ADDITIONAL INFORMATION

Product Categories/Family for NT5C3 antibody pair kit	Enzyme & Kinase
---	-------------------------------------

NCBI AND UNIPROT INFORMATION

NCBI GI #	51105879
NCBI GeneID	51251
Molecular Weight	32,559 Da
NCBI Official Full Name	5'-nucleotidase, cytosolic III
NCBI Official Synonym Full Names	5'-nucleotidase, cytosolic IIIA
NCBI Official Symbol	NT5C3A
NCBI Official Synonym Symbols	p36; PN-I; POMP; PSN1; UMPH; NT5C3; P5N-1; UMPH1; hUMP1; P5'N-1; cN-III
NCBI Protein Information	cytosolic 5'-nucleotidase 3A
UniProt Protein Name	Cytosolic 5'-nucleotidase 3A
UniProt Gene Name	NT5C3A
UniProt Synonym Gene Names	NT5C3; P5N1; UMPH1; 7-methylguanosine nucleotidase; cN-III; P5'N-1; P5N-1; PN-I

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/human-antibody-pair-kits/5-nucleotidase-cytosolic-iii-nt5c3/148907>