

## Human FDP ELISA Kit | FDP elisa kit

Human FDP (Fibrinogen Degradation Product) ELISA Kit

CATALOG #  
**AAA163736**

REACTIVITY  
**Human**


### PRODUCT SYNONYMS

FDP; Human FDP (Fibrinogen Degradation Product) ELISA Kit; FDP elisa kit

### PRODUCT SPECIFICATIONS

|                         |                                                                                                                                                                                                                                                                                                                                                                                            |
|-------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity              | Human                                                                                                                                                                                                                                                                                                                                                                                      |
| Specificity             | This assay has high sensitivity and excellent specificity for detection of Human FDP. No significant cross-reactivity or interference between Human FDP and analogues was observed. Note: Limited by current skills and knowledge, it is impossible for us to complete the cross- reactivity detection between Human FDP and all the analogues, therefore, cross reaction may still exist. |
| Sequence Length         | 128                                                                                                                                                                                                                                                                                                                                                                                        |
| Samples                 | Serum, Plasma, Biological Fluids                                                                                                                                                                                                                                                                                                                                                           |
| Assay Type              | Sandwich                                                                                                                                                                                                                                                                                                                                                                                   |
| Detection Range         | 1.563-100ng/mL                                                                                                                                                                                                                                                                                                                                                                             |
| Sensitivity             | Min: 0.938ng/mL; Max: 100ng/mL                                                                                                                                                                                                                                                                                                                                                             |
| Preparation and Storage | Store at 4 degree C.                                                                                                                                                                                                                                                                                                                                                                       |

### IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

|                                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|-----------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Related Product Information for FDP elisa kit | <p><b>Description:</b> The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of FDP in human serum, plasma and other biological fluids</p> <p><b>Principle of the Assay:</b> This ELISA kit uses Sandwich-ELISA as the method. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Human FDP. Standards or samples are added to the appropriate micro ELISA plate wells and bound by the specific antibody. Then a biotinylated detection antibody specific for Human FDP and Avidin-Horseradish Peroxidase (HRP) conjugate is added to each micro plate well successively and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Human FDP, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of a sulphuric acid solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm +/- 2 nm. The OD value is proportional to the concentration of Human FDP. You can calculate the concentration of Human FDP in the samples by comparing the OD of the samples to the standard curve.</p> |
|-----------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

NCBI AND UNIPROT INFORMATION

|                            |                           |
|----------------------------|---------------------------|
| NCBI GI #                  | 11991844                  |
| Molecular Weight           | 14,332 Da                 |
| NCBI Official Full Name    | fibrocyte-derived protein |
| UniProt Protein Name       | Otoraplin                 |
| UniProt Gene Name          | OTOR                      |
| UniProt Synonym Gene Names | FDP; MIAL                 |
| UniProt Entry Name         | OTOR_HUMAN                |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/fdp/163736>