

Canine TNF-alpha ELISA Kit | TNF-alpha elisa kit

Canine TNF-alpha (Tumor Necrosis Factor Alpha) ELISA Kit

CATALOG #
AAA165684

REACTIVITY
Canine


PRODUCT SYNONYMS

TNF-alpha; Canine TNF-alpha (Tumor Necrosis Factor Alpha) ELISA Kit; TNF-alpha elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Specificity	This kit recognizes Canine TNF-alpha in samples. No significant cross-reactivity or interference between Canine TNF-alpha and analogues was observed.
Sequence Length	233
Assay Type	Quantitative Sandwich
Samples	Serum, plasma and other biological fluids
Detection Range	15.625-1000pg/mL
Sensitivity	9.375pg/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, mid-range and high-level Canine TNF-alpha were tested 20 times on one plate, respectively.
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, mid-range and high-level Canine TNF-alpha were tested on 3 different plates, 20 replicates in each plate.
Preparation and Storage	Store at 4 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TNF-alpha elisa kit	Principle of the Assay: This ELISA kit uses the Sandwich-ELISA principle. The micro-ELISA plate provided in this kit has been pre-coated with an antibody specific to Canine TNF-alpha. Standards or samples are added to the micro-ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Canine TNF-alpha and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Canine TNF-alpha, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm ± 2 nm. The OD value is proportional to the concentration of Canine TNF-alpha. You can calculate the concentration of Canine TNF-alpha in the samples by comparing the OD of the samples to the standard curve. Intended Uses: This ELISA kit applies to the in vitro quantitative determination of Canine TNF-alpha concentrations in serum, plasma and other biological fluids.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	37210
NCBI GeneID	7124
UniProt Accession #	P01375; O43647; Q9P1Q2; Q9UIV3
Molecular Weight	25,644 Da
NCBI Official Full Name	TNF-alpha
NCBI Official Synonym Full Names	tumor necrosis factor
NCBI Official Symbol	TNF
NCBI Official Synonym Symbols	DIF; TNFA; TNFSF2; TNF-alpha
NCBI Protein Information	tumor necrosis factor; APC1 protein; TNF, macrophage-derived; TNF, monocyte-derived; TNF-a; cachectin; tumor necrosis factor ligand superfamily member 2; tumor necrosis factor-alpha
UniProt Protein Name	Tumor necrosis factor
UniProt Gene Name	TNF
UniProt Synonym Gene Names	TNFA; TNFSF2; TNF-a; NTF; ICD1; ICD2
UniProt Entry Name	TNFA_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/canine-elisa-kits/tnf-alpha/165684>