

## Human NSE (Neuron-Specific Enolase) CLIA Kit | NSE clia kit

Human NSE (Neuron-Specific Enolase) CLIA Kit

CATALOG #  
**AAA171217**

REACTIVITY  
**Human**


### PRODUCT SYNONYMS

NSE (Neuron-Specific Enolase); Human NSE (Neuron-Specific Enolase) CLIA Kit; Enolase gamma; Enolase g; Enolase 2; NSE clia kit

### PRODUCT SPECIFICATIONS

|                 |                                                                                                                                          |
|-----------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity      | Human                                                                                                                                    |
| Specificity     | This kit recognizes natural and recombinant NSE. No significant cross-reactivity or interference between NSE and analogues was observed. |
| Sequence Length | 434                                                                                                                                      |
| Assay Type      | Sandwich                                                                                                                                 |
| Detection Range | 0.469-30ng/mL                                                                                                                            |
| Sensitivity     | Min: 0.281ng/mL; Max: 30ng/mL                                                                                                            |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for NSE clia kit

**Intended Uses:** This CLIA kit applies to the in vitro quantitative determination of NSE concentrations in serum, plasma and other biological fluids.

**Principle of the Assay:** This kit uses Sandwich-CLIA as the method. The micro CLIA plate provided in this kit has been pre-coated with an antibody specific to NSE. Standards or samples are added to the appropriate micro CLIA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for NSE and Avidin-Horseradish Peroxidase (HRP) conjugate is added to each micro plate well successively and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain NSE, biotinylated detection antibody and Avidin-HRP conjugate will appear fluorescence. The Relative light unit (RLU) value is measured spectrophotometrically by the Chemiluminescence immunoassay analyzer. The RLU value is positively associated with the concentration of NSE. You can calculate the concentration of NSE in the samples by comparing the RLU value of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

|                               |                             |
|-------------------------------|-----------------------------|
| NCBI GI #                     | 359719155                   |
| NCBI GeneID                   | 102190860                   |
| Molecular Weight              | 47,266 Da                   |
| NCBI Official Full Name       | neuron-specific enolase     |
| NCBI Official Symbol          | ENO2                        |
| NCBI Official Synonym Symbols | NSE                         |
| NCBI Protein Information      | enolase 2 (gamma, neuronal) |
| UniProt Protein Name          | Neuron-specific enolase     |
| UniProt Gene Name             | NSE                         |
| UniProt Entry Name            | G9JNA8_CAPHI                |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://www.aaabiotech.com/human-clia-kits/nse-neuron-specific-enolase/171217>