

## 2-Step Plus Poly-HRP Detection Kit

2-Step Plus Poly-HRP Anti Rabbit/Mouse IgG Detection System (With DAB Solution)

CATALOG #  
**AAA174691**

APPLICATIONS  
**Immunohistochemistry**

### PRODUCT SYNONYMS

2-Step Plus Poly-HRP; 2-Step Plus Poly-HRP Anti Rabbit/Mouse IgG Detection System (With DAB Solution); 2-Step Plus Poly-HRP detection kit

### PRODUCT SPECIFICATIONS

|                                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|----------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Applicable Applications for 2-Step Plus Poly-HRP detection kit | IHC (Immunohistochemistry)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Sample dyeing                                                  | <ol style="list-style-type: none"> <li>1. Dewax and hydrate the paraffin section.</li> <li>2. Repair antigen of the tissue section according to special requirements of applied primary antibody.</li> <li>3. Incubate with (3% H<sub>2</sub>O<sub>2</sub>) for 10min to eliminate endogenous peroxidase activity.</li> <li>4. Wash with PBS or TBS, 2 min × 3 times.</li> <li>5. Drain the PBS with absorbent paper, then add (Normal Goat Serum) to the section. Incubate for 30 min at 37°C.</li> <li>6. Dry the liquid around the section with absorbent paper, and draw a circle around the tissue with an oily pen. Add primary antibody with proper dilution ratio, incubate at room temperature or 37°C for 1~2h or at 4°C overnight (then rewarm at 37°C for 30 min). Wash with PBS or TBS, 2 min × 3 times.</li> <li>7. Add (Polymer Helper), incubate at room temperature or 37°C for 20min. Wash with PBS or TBS, 2 min × 3 times.</li> <li>8. Add (Polyperoxidase-anti-Mouse/Rabbit IgG), incubate at room temperature or 37°C for 20~30min. Wash with PBS or TBS, 2 min × 3 times.</li> <li>9. Add 1 drop (approximately 50 uL) of (DAB Concentrate) into each 1 mL of (DAB Substrate), mix fully and the mixed reagent is the DAB Working Solution. Prepare fresh solution before use and the prepared solution should be stored in the dark. Fresh prepared DAB Working Solution is valid within 4 hours and the unused solution must be abandoned.</li> <li>10. Take control of the DAB coloration period, the color of tan or brownish yellow is the positive signal. Avoid of excessive reaction.</li> <li>11. Wash the section with deionized water terminate the chromogenic reaction, then operate the procedures of counterstaining, dehydrating, transparentizing and sealing.</li> </ol> |
| Preparation and Storage                                        | Store at 2~8°C, shading light. Avoid of freezing. Valid for 12 months. The reagents are valid within 6 months after opening.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |

### ADDITIONAL INFORMATION

|                                                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|--------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Related Product Information for 2-Step Plus Poly-HRP detection kit | <p>Background: 2-step plus is a two-step immunohistochemical broad spectrum detection reagent. It polymerizes monovalent Fab fragments of secondary antibody and enzyme, which replaces the secondary antibody and tertiary antibody in traditional method, can directly amplify the binding signal of antibody-antigen. This method not only retains the specific binding ability of antibody with antigen, but also can effectively avoid space steric hindrance caused by excessive polymer molecules. Compared with the traditional SP three-step method, this kit has the characteristics of simple, rapid, high-sensitivity. This system abandons the using of biotin, so it can avoid background staining by endogenous biotin. It can be used in IHC, in which the primary antibody is monoclonal/polyclonal antibody derived from rabbit/mouse. The unique polymer auxiliary agent for macromolecular detection system can detect the binding primary antibody better. Diluent for DAB concentrated solution has been provided in this kit to avoid the influence of the different water acidity and alkalinity on the DAB Chromogenic Agent.</p> |
|--------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.