

## Human Sonic hedgehog (SHH) ELISA Kit | SHH elisa kit

Human Sonic hedgehog (SHH) ELISA Kit

CATALOG #  
**AAA180107**

REACTIVITY  
**Human**

### PRODUCT SYNONYMS

Sonic hedgehog (SHH); Human Sonic hedgehog (SHH) ELISA Kit; SHH elisa kit

### PRODUCT SPECIFICATIONS

|                         |                                       |
|-------------------------|---------------------------------------|
| Reactivity              | Human                                 |
| Specificity             | No cross-reaction with other factors. |
| Sequence Length         | 418                                   |
| Assay Type              | Sandwich                              |
| Detection Range         | 10 ng/mL-0.156 ng/mL                  |
| Sensitivity             | Up to 0.05 ng/mL                      |
| Intra-assay Precision   | <= 8%                                 |
| Inter-assay Precision   | <= 12%                                |
| Preparation and Storage | Store all reagents at 2-8 degree C.   |

### ADDITIONAL INFORMATION

#### Related Product Information for SHH elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Human SHH monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

### NCBI AND UNIPROT INFORMATION

|                                  |                                            |
|----------------------------------|--------------------------------------------|
| NCBI GI #                        | <a href="#">1171142</a>                    |
| NCBI GeneID                      | <a href="#">30269</a>                      |
| Molecular Weight                 | 46,403 Da                                  |
| NCBI Official Full Name          | sonic hedgehog                             |
| NCBI Official Synonym Full Names | sonic hedgehog a                           |
| NCBI Official Symbol             | shha                                       |
| NCBI Official Synonym Symbols    | shh; syu; vhh1; vhh-1; fc83d08; wu:fc83d08 |
| NCBI Protein Information         | sonic hedgehog protein A; sonic-you        |
| UniProt Protein Name             | Sonic hedgehog protein A                   |

|                            |                 |
|----------------------------|-----------------|
| UniProt Gene Name          | shha            |
| UniProt Synonym Gene Names | shh; vhh1; SHHA |
| UniProt Entry Name         | SHH_DANRE       |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/sonic-hedgehog-shh/180107>