

Human NOGO-A ELISA Kit | NOGO-A elisa kit

Human NOGO-A ELISA Kit

CATALOG #
AAA180891

REACTIVITY
Human


PRODUCT SYNONYMS

NOGO-A; Human NOGO-A ELISA Kit; NOGO-A elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	294
Assay Type	Sandwich
Detection Range	20 ng/mL-0.312 ng/mL
Sensitivity	Up to 0.06 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for NOGO-A elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Human NOGO-A monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	76879695
NCBI GeneID	57142
Molecular Weight	108,450 Da

NCBI Official Full Name	NOGO-A, partial
NCBI Official Synonym Full Names	reticulon 4
NCBI Official Symbol	RTN4
NCBI Official Synonym Symbols	ASY; NSP; NOGO; NOGOC; RTN-X; NOGO-A; NSP-CL; Nogo-B; Nogo-C; RTN4-A; RTN4-C; RTN4-B1; RTN4-B2; NI220/250; Nbla00271; Nbla10545
NCBI Protein Information	reticulon-4
UniProt Protein Name	Reticulon-4
UniProt Gene Name	RTN4
UniProt Synonym Gene Names	KIAA0886; NOGO; Nogo protein; NSP
UniProt Entry Name	RTN4_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/nogo-a/180891>