

Canine S100 Calcium Binding Protein ELISA Kit | S100 elisa kit

Canine S100 Calcium Binding Protein (S100) ELISA Kit

CATALOG #
AAA181043

REACTIVITY
Canine


PRODUCT SYNONYMS

S100 Calcium Binding Protein; Canine S100 Calcium Binding Protein (S100) ELISA Kit; S100 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Specificity	No cross-reaction with other factors.
Sequence Length	92
Samples	Canine Serum, Plasma Or Cell Culture Supernatant And Organizations In The Natural And Recombinant Rln Concentration
Assay Type	Sandwich
Detection Range	20 ng/ml-0.312 ng/ml
Sensitivity	12 pg/ml.
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for S100 elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Canine RLN monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #

12804681

NCBI GeneID	6285
Molecular Weight	10,713 Da
NCBI Official Full Name	S100 calcium binding protein B
NCBI Official Synonym Full Names	S100 calcium binding protein B
NCBI Official Symbol	S100B
NCBI Official Synonym Symbols	NEF; S100; S100-B; S100beta
NCBI Protein Information	protein S100-B
UniProt Protein Name	Protein S100-B
UniProt Gene Name	S100B
UniProt Entry Name	S100B_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/canine-elisa-kits/s100-calcium-binding-protein/181043>