

Bovine Desmin (Des) ELISA Kit | Des elisa kit

Bovine Desmin(Des) ELISA Kit

CATALOG # AAA181267	REACTIVITY Bovine
-------------------------------	-----------------------------

PRODUCT SYNONYMS

Desmin (Des); Bovine Desmin(Des) ELISA Kit; Des elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Bovine
Sequence Length	470
Assay Type	Sandwich
Detection Range	200 ng/mL-3.12 ng/mL
Sensitivity	Up to 0.6 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

ADDITIONAL INFORMATION

Related Product Information for Des elisa kit	Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Bovine Des monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	55749932
NCBI GeneID	1674
NCBI Accession #	NP_001918.3
NCBI GenBank Nucleotide #	NM_001927.3
Molecular Weight	53,536 Da
NCBI Official Full Name	desmin
NCBI Official Synonym Full Names	desmin
NCBI Official Symbol	DES
NCBI Official Synonym Symbols	CSM1; CSM2; LGMD2R
NCBI Protein Information	desmin

UniProt Protein Name	Desmin
UniProt Gene Name	DES
UniProt Entry Name	DESM_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/bovine-ovine-elisa-kits/desmin-des/181267>