

## Human Kallikrein 12 (KLK12) ELISA Kit | KLK12 elisa kit

Human Kallikrein 12 (KLK12) ELISA Kit

CATALOG #  
AAA182370

REACTIVITY  
Human


### PRODUCT SYNONYMS

Kallikrein 12 (KLK12); Human Kallikrein 12 (KLK12) ELISA Kit; KLK-L5; KLKL5; Kallikrein-Related Peptidase 12; Kallikrein-like protein 5; KLK12 elisa kit

### PRODUCT SPECIFICATIONS

|                                                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|---------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity                                        | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Specificity                                       | This assay has high sensitivity and excellent specificity for detection of Kallikrein 12 (KLK12). No significant cross-reactivity or interference between Kallikrein 12 (KLK12) and analogues was observed.                                                                                                                                                                                                                                                                                |
| Sequence Length                                   | 111                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Assay Type                                        | Competitive Inhibition                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Samples                                           | Serum, Plasma, Tissue Homogenates and Other Biological Fluids                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Detection Range                                   | 12.35-1000ng/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Sensitivity                                       | < 7.14ng/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Application                                       | Enzyme-linked immunosorbent assay for Antigen Detection.                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Intra-assay Precision (Precision within an assay) | 3 samples with low, middle and high level Kallikrein 12 (KLK12) were tested 20 times on one plate, respectively.                                                                                                                                                                                                                                                                                                                                                                           |
| Inter-assay Precision (Precision between assays)  | 3 samples with low, middle and high level Kallikrein 12 (KLK12) were tested on 3 different plates, 8 replicates in each plate.                                                                                                                                                                                                                                                                                                                                                             |
| CV(%) =                                           | SD/meanX100                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Intra-Assay                                       | CV<10%                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Inter-Assay                                       | CV<12%                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Preparation and Storage                           | The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.<br>To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end. |

### IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

|                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|-------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Related Product Information for KLK12 elisa kit | This assay employs the competitive inhibition enzyme immunoassay technique. A monoclonal antibody specific to Kallikrein 12 (KLK12) has been pre-coated onto a microplate. A competitive inhibition reaction is launched between biotin labeled Kallikrein 12 (KLK12) and unlabeled Kallikrein 12 (KLK12) (Standards or samples) with the pre-coated antibody specific to Kallikrein 12 (KLK12). After incubation the unbound conjugate is washed off. Next, avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. The amount of bound HRP conjugate is reverse proportional to the concentration of Kallikrein 12 (KLK12) in the sample. After addition of the substrate solution, the intensity of color developed is reverse proportional to the concentration of Kallikrein 12 (KLK12) in the sample. |
|-------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

NCBI AND UNIPROT INFORMATION

|                                  |                                 |
|----------------------------------|---------------------------------|
| NCBI GI #                        | 22208989                        |
| NCBI GeneID                      | 43849                           |
| NCBI Accession #                 | NP_665902.1                     |
| NCBI GenBank Nucleotide #        | NM_145895.1                     |
| UniProt Accession #              | Q9UKR0; Q9UKR1; Q9UKR2          |
| Molecular Weight                 | 11,984 Da                       |
| NCBI Official Full Name          | kallikrein-12 isoform 3         |
| NCBI Official Synonym Full Names | kallikrein related peptidase 12 |
| NCBI Official Symbol             | KLK12                           |
| NCBI Official Synonym Symbols    | KLKL5; KLK-L5                   |
| NCBI Protein Information         | kallikrein-12                   |
| UniProt Protein Name             | Kallikrein-12                   |
| UniProt Gene Name                | KLK12                           |
| UniProt Synonym Gene Names       | KLKL5; KLK-L5                   |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.