

Nucleoprotein (N) Recombinant Protein | N recombinant protein

Recombinant Crimean-Congo hemorrhagic fever virus Nucleoprotein (N)

CATALOG # AAA18641	HOST E Coli or Yeast or Baculovirus or Mammalian Cell
------------------------------	---


PRODUCT SYNONYMS

Nucleoprotein (N); Recombinant Crimean-Congo hemorrhagic fever virus Nucleoprotein (N); Nucleoprotein; Nucleocapsid protein; Protein N; N recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	1-482aa; Full Length
Sequence	MENKIEVNNK DEMNKFEEF KKGNGLVDTF TNPYSFCESV PNLERFVFQM ASATDDAQKD SIYASALVEA TKFCAPIYEC AWWSTGIVK KGLEWFEKNA GTIKSWDESY IELKVEVPKI EQLANYQQAA LKWRKDIGFR VNANTAALSH KVLAEYKVPK EIVMSVKEML SDMIRRRNLI LNRGGDENPR GPVSREHVEW CREFVKGKYI MAFNPPWGD I NKSGRSGIAL VATGLAKLAE TEGKGVFDEA KKTVEALNGY LDKHKDEVK ASADNMITNL LKHIKAQEL YKNSSALRAQ GAQIDTAFSS YYWLYKAGVT PETFPTVSQF LFELGKQPRG TKKMKKALLS TPMKWGKKLY ELFADDSFQQ NRIYMHPAVL TAGRISEMGV CFGTIPVANP DDAAQGSQHT KSILNLRNT ETNNPCAkti VKLFEIQKTG FNIQDMDIVA SEHLLHQSLV GKQSPFQNAV NVKGNATSAN II
Species	Crimean-Congo hemorrhagic fever virus (isolate C68031) (CCHFV)
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

IMAGES



SDS-PAGE

ADDITIONAL INFORMATION

Related Product Information for N recombinant protein	Recombinant Crimean-Congo hemorrhagic fever virus Nucleoprotein (N)
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	127869
UniProt Accession #	P27317
Molecular Weight	68 kDa
NCBI Official Full Name	Nucleoprotein
UniProt Protein Name	Nucleoprotein
UniProt Gene Name	N
UniProt Synonym Gene Names	Protein N
UniProt Entry Name	NCAP_CCHFV

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/n-recombinant-protein-e-coli/18641>