

High affinity transport system protein p37 (p37) Recombinant Protein | p37 recombinant protein

Recombinant Mycoplasma hyorhinis High affinity transport system protein p37 (p37)

CATALOG # AAA18655	HOST E Coli or Yeast or Baculovirus or Mammalian Cell
------------------------------	---


PRODUCT SYNONYMS

High affinity transport system protein p37 (p37); Recombinant Mycoplasma hyorhinis High affinity transport system protein p37 (p37); Mycoplasma hyorhinis; p37 recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	24-403aa; Full Length of Mature Protein
Sequence	CSNTGVVKQE DVSVSQGQWD KSITFGVSEA WLNKKKGGEK VNKEVINTFL ENFKKEFNKL KNANDKTKNF DDVDFKVTP I QDFTVLLNNL STDNPELDFG INASGKLVEF LKNNPGIITP ALETTTNSFV FDKKDKFYV DGTDSPLVK IAKEINKIFV ETPYASWTDE NHKWNGNVYQ SVYDPTVQAN FYRGMWIKG NDETLAKIKK AWNDKDWNTF RNFGILHGKD NSFSKFKLEE TILKNHFQNK FTTLNEDRSA HPNAYKQKSA DTLGLDDFH IAFSEEGSFA WTHNKSATKP FETKANEME ALIVTNPIPY DVGVRKSVN QLEQNLIVQT FINLAKNKQD TYGPLLGYN G YKKIDNFQKE IVEVEKAIAK
Species	Mycoplasma hyorhinis
Preparation and Storage	Store at -20 degree C, for extended storage, conserve at -20 degree C or -80 degree C.

IMAGES



SDS-PAGE

ADDITIONAL INFORMATION

Related Product Information for p37 recombinant protein	P37 is part of a high-affinity transport system.
---	--

References	"A mycoplasma high-affinity transport system and the in vitro invasiveness of mouse sarcoma cells." Dudler R., Schmidhauser C., Parish R.W., Wettenhall R.E.H., Schmidt T. EMBO J. 7:3963-3970(1988)
------------	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	129354
UniProt Accession #	P15363
Molecular Weight	59.5 kDa
NCBI Official Full Name	High affinity transport system protein p37
UniProt Protein Name	High affinity transport system protein p37
UniProt Gene Name	p37
UniProt Entry Name	P37_MYCHR

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
 Online datasheet: <https://www.aaabiotech.com/p37-recombinant-protein-e-coli/18655>