

Human Vitamin K Epoxide Reductase Complex Subunit 1 (VKORC1) ELISA Kit | VKORC1 elisa kit

Human Vitamin K Epoxide Reductase Complex Subunit 1 (VKORC1) ELISA Kit

CATALOG # AAA187192	REACTIVITY Human
-------------------------------	----------------------------


PRODUCT SYNONYMS

Vitamin K Epoxide Reductase Complex Subunit 1 (VKORC1); Human Vitamin K Epoxide Reductase Complex Subunit 1 (VKORC1) ELISA Kit; VKCFD2; VKOR; Vitamin K Dependent Clotting Factors Deficiency 2; Vitamin K1 2; 3-epoxide reductase subunit 1; VKORC1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Vitamin K Epoxide Reductase Complex Subunit 1 (VKORC1). No significant cross-reactivity or interference between Vitamin K Epoxide Reductase Complex Subunit 1 (VKORC1) and analogues was observed.
Sequence Length	92
Assay Type	Double-antibody Sandwich
Samples	Tissue Homogenates, Cell Lysates and Other Biological Fluids
Detection Range	1.56-100ng/mL
Sensitivity	< 0.60ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Vitamin K Epoxide Reductase Complex Subunit 1 (VKORC1) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Vitamin K Epoxide Reductase Complex Subunit 1 (VKORC1) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for VKORC1 elisa kit

The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Vitamin K Epoxide Reductase Complex Subunit 1 (VKORC1). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Vitamin K Epoxide Reductase Complex Subunit 1 (VKORC1). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Vitamin K Epoxide Reductase Complex Subunit 1 (VKORC1), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Vitamin K Epoxide Reductase Complex Subunit 1 (VKORC1) in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	62511226
NCBI GeneID	79001
UniProt Accession #	Q9BQB6; Q6UX90; Q7Z2R4; A6NIQ6; B2R4Z6
Molecular Weight	9,875 Da
NCBI Official Full Name	Vitamin K epoxide reductase complex subunit 1
NCBI Official Synonym Full Names	vitamin K epoxide reductase complex subunit 1
NCBI Official Symbol	VKORC1
NCBI Official Synonym Symbols	VKOR; MST134; MST576; VKCFD2; EDTP308
NCBI Protein Information	vitamin K epoxide reductase complex subunit 1
UniProt Protein Name	Vitamin K epoxide reductase complex subunit 1
UniProt Gene Name	VKORC1
UniProt Synonym Gene Names	VKOR

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/vitamin-k-epoxide-reductase-complex-subunit-1-vkorc1/187192>