

Human Trans Golgi Network Protein 2 (TGOLN2) ELISA Kit | TGOLN2 elisa kit

Human Trans Golgi Network Protein 2 (TGOLN2) ELISA Kit

CATALOG #
AAA187230

REACTIVITY
Human


PRODUCT SYNONYMS

Trans Golgi Network Protein 2 (TGOLN2); Human Trans Golgi Network Protein 2 (TGOLN2) ELISA Kit; TGN51; TGN46; TGN48; TGN38; TTN2; Trans-Golgi Network Protein; TGOLN2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Trans Golgi Network Protein 2 (TGOLN2). No significant cross-reactivity or interference between Trans Golgi Network Protein 2 (TGOLN2) and analogues was observed.
Sequence Length	437
Assay Type	Double-antibody Sandwich
Samples	Tissue Homogenates, Cell Lysates and Other Biological Fluids
Detection Range	0.312-20ng/mL
Sensitivity	< 0.112ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Trans Golgi Network Protein 2 (TGOLN2) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Trans Golgi Network Protein 2 (TGOLN2) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TGOLN2 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Trans Golgi Network Protein 2 (TGOLN2). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Trans Golgi Network Protein 2 (TGOLN2). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Trans Golgi Network Protein 2 (TGOLN2), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Trans Golgi Network Protein 2 (TGOLN2) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	14250103
NCBI GeneID	10618
Molecular Weight	29,541 Da
NCBI Official Full Name	Trans-golgi network protein 2
NCBI Official Synonym Full Names	trans-golgi network protein 2
NCBI Official Symbol	TGOLN2
NCBI Official Synonym Symbols	TGN38; TGN46; TGN48; TGN51; TTGN2
NCBI Protein Information	trans-Golgi network integral membrane protein 2
UniProt Protein Name	Trans-Golgi network integral membrane protein 2
UniProt Gene Name	TGOLN2
UniProt Synonym Gene Names	TGN46; TGN51

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/trans-golgi-network-protein-2-tgoln2/187230>