

## Rat Myoglobin (MYO) CLIA Kit | MYO clia kit

Rat Myoglobin (MYO) CLIA Kit

|                               |                          |                                                             |
|-------------------------------|--------------------------|-------------------------------------------------------------|
| CATALOG #<br><b>AAA188918</b> | REACTIVITY<br><b>Rat</b> | APPLICATIONS<br><b>CLIA (Chemiluminescence Immunoassay)</b> |
|-------------------------------|--------------------------|-------------------------------------------------------------|


### PRODUCT SYNONYMS

Myoglobin (MYO); Rat Myoglobin (MYO) CLIA Kit; Myoglobin; PVALB; MB; MYO clia kit

### PRODUCT SPECIFICATIONS

|                                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity                               | Rat                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Specificity                              | This assay has high sensitivity and excellent specificity for detection of Myoglobin (MYO). No significant cross-reactivity or interference between Myoglobin (MYO) and analogues was observed.                                                                                                                                                                                                                                                                                                                                                                                  |
| Sequence Length                          | 154                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Assay Type                               | Double-antibody Sandwich                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Samples                                  | Serum, plasma, tissue homogenates and other biological fluids                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Detection Range                          | 0.14-100ng/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Sensitivity                              | The minimum detectable dose of this kit is typically less than 0.05ng/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Applicable Applications for MYO clia kit | CLIA (Chemiluminescence Immunoassay)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Precision                                | Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level Myoglobin (MYO) were tested 20 times on one plate, respectively.<br>Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level Myoglobin (MYO) were tested on 3 different plates, 8 replicates in each plate.<br>CV(%) = SD/meanX100<br>Intra-Assay: CV<10%<br>Inter-Assay: CV<12%                                                                                                                                                                 |
| Test Principle                           | The microplate provided in this kit has been pre-coated with an antibody specific to Myoglobin (MYO). Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to Myoglobin (MYO). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then the mixture of substrate A and B is added to generate glow light emission kinetics. Upon plate development, the intensity of the emitted light is proportional to the Myoglobin (MYO) level in the sample or standard. |
| Preparation and Storage                  | The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.<br>To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.                                                                                       |

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

|                           |             |
|---------------------------|-------------|
| NCBI GI #                 | 11024650    |
| NCBI Accession #          | NP_067599.1 |
| NCBI GenBank Nucleotide # | NP_067599.1 |
| UniProt Accession #       | Q9QZ76      |
| NCBI Official Full Name   | myoglobin   |
| UniProt Protein Name      | Myoglobin   |
| UniProt Gene Name         | Mb          |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://www.aaabiotech.com/rat-clia-kits/myoglobin-myo/188918>