

Mouse Estrone (E1) ELISA Kit | E1 elisa kit

Mouse Estrone (E1) ELISA Kit

CATALOG #
AAA190328

REACTIVITY
Mouse

PRODUCT SYNONYMS

Estrone (E1); Mouse Estrone (E1) ELISA Kit; E1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of Mouse E1. No significant cross-reactivity or interference between Mouse E1 and analogues was observed.
Assay Type	Quantitative Competitive
Samples	Serum, plasma and other biological fluids
Detection Range	80 pg/mL - 2000 pg/mL
Sensitivity	40 pg/mL (mean of 6 independent assays)
Intra-assay Precision	Intra-assay Precision (Precision within an assay) Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision. Intra-Assay: CV<8%
Inter-assay Precision	Inter-assay Precision (Precision between assays) Three samples of known concentration were tested in forty separate assays to assess inter-assay precision. CV(%) = SD/meanX100. Inter-Assay: CV<10%

ADDITIONAL INFORMATION

Related Product Information for E1 elisa kit

Principle of the Assay: This assay employs the competitive enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with goat-anti-rabbit antibody. Standards or samples are added to the appropriate microtiter plate wells with an antibody specific for E1 and Horseradish Peroxidase (HRP) conjugated E1. The competitive inhibition reaction is launched between with HRP labeled E1 and unlabeled E1 with the antibody. A substrate solution is added to the wells and the color develops in opposite to the amount of E1 in the sample. The color development is stopped and the intensity of the color is measured.

Intended Uses: For the quantitative detection of Mouse Estrone (E1) concentration in serum, plasma and other biological fluids.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/mouse-elisa-kits/estrone-e1/190328>