

Human D-Aspartate Oxidase (DDO) ELISA Kit | DDO elisa kit

Human D-Aspartate Oxidase (DDO) ELISA Kit

CATALOG #
AAA202084

REACTIVITY
Human


PRODUCT SYNONYMS

D-Aspartate Oxidase (DDO); Human D-Aspartate Oxidase (DDO) ELISA Kit; DDO elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	341
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	1ng/ml-16ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for DDO elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of DDO in the sample, This DDO ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus DDO concentration. The concentration of DDO in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	261278134
NCBI GeneID	8528
NCBI Official Full Name	D-aspartate oxidase

NCBI Official Synonym Full Names	D-aspartate oxidase
NCBI Official Symbol	DDO
NCBI Official Synonym Symbols	DASOX; DDO-1; DDO-2
NCBI Protein Information	D-aspartate oxidase

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/d-aspartate-oxidase-ddo/202084>