

## Human Tartrate-Resistant Acid Phosphatase (TRAP) ELISA Kit | TRAP elisa kit

Human Tartrate-Resistant Acid Phosphatase (TRAP) ELISA Kit

CATALOG #  
**AAA202748**

REACTIVITY  
**Human**

### PRODUCT SYNONYMS


Tartrate-Resistant Acid Phosphatase (TRAP); Human Tartrate-Resistant Acid Phosphatase (TRAP) ELISA Kit; TRAP elisa kit

### PRODUCT SPECIFICATIONS

Reactivity

Human

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for TRAP elisa kit

**Intended Uses:** This TRAP ELISA kit is intended Laboratory for Research use only. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of TRAP in the sample, this TRAP ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus TRAP concentration. The concentration of TRAP in the samples is then determined by comparing the O.D. of the samples to the standard curve.

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)

Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/tartrate-resistant-acid-phosphatase-trap/202748>