

Human VE-Cadherin ELISA Kit | CDH5 elisa kit

Human VE-Cadherin ELISA Kit

CATALOG #
AAA203364

REACTIVITY
Human

PRODUCT SYNONYMS

VE-Cadherin; Human VE-Cadherin ELISA Kit; CDH5 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	784
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	0.5 mug/ml-8 mug/ml
Sensitivity	0.1 mug/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CDH5 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of VE-CADHERIN in the sample, This VE-CADHERIN ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus VE-CADHERIN concentration. The concentration of VE-CADHERIN in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	166362713
NCBI GeneID	1003

NCBI Accession #	NP_001786.2
NCBI GenBank Nucleotide #	NP_001786.2
UniProt Accession #	P33151
NCBI Official Full Name	cadherin-5 preproprotein
NCBI Official Synonym Full Names	cadherin 5
NCBI Official Symbol	CDH5
NCBI Official Synonym Symbols	7B4; CD144
NCBI Protein Information	cadherin-5
UniProt Protein Name	Cadherin-5
UniProt Gene Name	CDH5
UniProt Synonym Gene Names	VE-cadherin
UniProt Entry Name	CADH5_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/ve-cadherin/203364>