

Human Japanese Encephalitis IgM ELISA Kit

Human Japanese Encephalitis IgM Antibody ELISA Kit

CATALOG #
AAA203696

REACTIVITY
Human

PRODUCT SYNONYMS

Japanese Encephalitis IgM; Human Japanese Encephalitis IgM Antibody ELISA Kit; Japanese Encephalitis IgM elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Intra-assay Precision	Intra-assay CV is less than 15%.
Inter-assay Precision	Inter-assay CV is less than 15%.

ADDITIONAL INFORMATION

Related Product Information for Japanese Encephalitis IgM elisa kit

Intended Uses: This Japanese encephalitis IgM antibody ELISA kit is intended Laboratory for Research use only!. In order to determine whether contains Japanese encephalitis IgM antibody in the sample, this Japanese encephalitis IgM antibody ELISA Kit includes Negative control and Positive control. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. The color depth was positively correlated with the Japanese encephalitis IgM antibody in the sample.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/japanese-encephalitis-igm/203696>