

## Human Specific IgE4 Against Dust Mite ELISA Kit

Human Specific IgE4 Against Dust Mite Antibody ELISA Kit

CATALOG #  
**AAA203787**

REACTIVITY  
**Human**

### PRODUCT SYNONYMS

Specific IgE4 Against Dust Mite; Human Specific IgE4 Against Dust Mite Antibody ELISA Kit; Specific IgE4 Against Dust Mite elisa kit

### PRODUCT SPECIFICATIONS

|                 |                                                                   |
|-----------------|-------------------------------------------------------------------|
| Reactivity      | Human                                                             |
| Samples         | Serum, plasma, cell culture supernate and other biological fluids |
| Assay Type      | Sandwich                                                          |
| Detection Range | 15 IU/mL-240 IU/mL                                                |
| Sensitivity     | 1 IU/ml                                                           |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for Specific IgE4 Against Dust Mite elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of SPECIFIC IGE4 AGAINST DUST MITE ANTIBODY in the sample, This SPECIFIC IGE4 AGAINST DUST MITE ANTIBODY ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus SPECIFIC IGE4 AGAINST DUST MITE ANTIBODY concentration. The concentration of SPECIFIC IGE4 AGAINST DUST MITE ANTIBODY in the samples is then determined by comparing the O.D. of the samples to the standard curve.

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.