

Human Varicella-Zoster Virus ELISA Kit

Human Varicella-Zoster Virus ELISA Kit

CATALOG #
AAA203800

REACTIVITY
Human

PRODUCT SYNONYMS

Varicella-Zoster Virus; Human Varicella-Zoster Virus ELISA Kit; Varicella-Zoster Virus elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Intra-assay Precision	Intra-assay CV is less than 15%.
Inter-assay Precision	Inter-assay CV is less than 15%.

ADDITIONAL INFORMATION

Related Product Information for Varicella-Zoster Virus elisa kit

Intended Uses: This Varicella-zoster virus ELISA kit is intended Laboratory for Research use only and is not for use in test or therapeutic procedures. In order to determine whether contains Varicella-zoster virus in the sample, this Varicella-zoster virus ELISA Kit includes Negative control and Positive control. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. The color depth was positively correlated with the Varicella-zoster virus in the sample.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/varicella-zoster-virus/203800>