

Human N. Gonorrhoeae ELISA Kit

Human N. Gonorrhoeae ELISA Kit

CATALOG #
AAA203986

REACTIVITY
Human

PRODUCT SYNONYMS

N. Gonorrhoeae; Human N. Gonorrhoeae ELISA Kit; N. Gonorrhoeae elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Intra-assay Precision	Intra-assay CV is less than 15%.
Inter-assay Precision	Inter-assay CV is less than 15%.

ADDITIONAL INFORMATION

Related Product Information for N. Gonorrhoeae elisa kit

Intended Uses: This N. gonorrhoeae ELISA kit is intended Laboratory for Research use only. In order to determine whether contains N. gonorrhoeae in the sample, this N. gonorrhoeae ELISA Kit includes Negative control and Positive control. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. The color depth was positively correlated with the N. gonorrhoeae in the sample.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/n-gonorrhoeae/203986>