

Human Total Procollagen Type 1 Intact N-Terminal Propeptide ELISA Kit

Human Total Procollagen Type 1 Intact N-Terminal Propeptide ELISA Kit

CATALOG #
AAA204547

REACTIVITY
Human


PRODUCT SYNONYMS

Total Procollagen Type 1 Intact N-Terminal Propeptide; Human Total Procollagen Type 1 Intact N-Terminal Propeptide ELISA Kit; Total Procollagen Type 1 Intact N-Terminal Propeptide elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	0.625ng/ml-10ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Total Procollagen Type 1 Intact N-Terminal Propeptide elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of TOTAL PROCOLLAGEN TYPE1 INTACTN-TERMINAL PROPEPTIDE in the sample, This TOTAL PROCOLLAGEN TYPE1 INTACTN-TERMINAL PROPEPTIDE ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus TOTAL PROCOLLAGEN TYPE1 INTACTN-TERMINAL PROPEPTIDE concentration. The concentration of TOTAL PROCOLLAGEN TYPE1 INTACTN-TERMINAL PROPEPTIDE in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.