

Human Tumor Necrosis Factor-alpha Converting Enzyme ELISA Kit

Human Tumor Necrosis Factor-alpha Converting Enzyme ELISA Kit

CATALOG #
AAA204601

REACTIVITY
Human


PRODUCT SYNONYMS

Tumor Necrosis Factor-alpha Converting Enzyme; Human Tumor Necrosis Factor-alpha Converting Enzyme ELISA Kit; Tumor Necrosis Factor-alpha Converting Enzyme elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	12.5 mU/ml-200 mU/ml
Sensitivity	1 mU/mL

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Tumor Necrosis Factor-alpha Converting Enzyme elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of TUMOR NECROSIS FACTOR-alpha CONVERTING ENZYME in the sample, This TUMOR NECROSIS FACTOR-alpha CONVERTING ENZYME ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus TUMOR NECROSIS FACTOR-alpha CONVERTING ENZYME concentration. The concentration of TUMOR NECROSIS FACTOR-alpha CONVERTING ENZYME in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.